XXL-Rehab Patientstol

XXL-Rehab Patient/Dinner chair

DK

UK

MAX 325 KG

DK

Kan med fordel anvendes til:

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Instructions for XXL-Rehab Patient/Dinner chair

⚠ CAUTION ⚠ CAUTION ⚠

- This product should not be used without proper instruction from a healthcare professional.
- Using the product for individuals weighing more than the maximum weight-bearing capacity may pose a risk for damaging the product.
- XXL-Rehab assumes no responsibility for any damage or injury caused by improper installation, assembly or use of this product.
- If components are damaged or missing, contact XXL-Rehab immediately. Do NOT use substitute parts.
- XXL-Rehab cannot be held liable for faults or accidents occurring after repairs by service staff without written authorisation from XXL-Rehab.

⚠ To ensure proper use of the XXL-Rehab Patient/Dinner chair, please read these instructions thoroughly before using the product. Pay particular attention to the warning notes marked with ⚠. Also remember that in order to prevent the risk of injury to the user, the patient/dinner chair should be adjusted only by trained caregivers or aides or by individuals who have been instructed in its proper use.

USE

The XXL-Rehab Patient/Dinner chair is designed to help an individual weighing up to 325 kg when seating is needed e.g. in a hospital, at a dinner table.

⚠ Using the product for individuals weighing more than 325 kg may pose a risk of tipping or damaging the patient/dinner chair.

⚠ Always use the XXL-Rehab Patient/Dinner chair on level, hard flooring with a maximum slope of 3 degrees, as otherwise there is a risk of tipping.

BEFORE USE

The bolts for adjusting the legs and backrest must be properly mounted and tightened in the holes on the side of the patient/dinner chair.

⚠ Do not use the patient/dinner chair before the bolts for adjusting the legs and backrest have been properly mounted and tightened, as otherwise there is a risk of the user falling off to the side of the chair.

DURING USE

⚠ To prevent the risk of tipping, only use the patient/dinner chair on a level, hard surface.

⚠ If the user requires assistance during transfers to the patient/dinner chair, the aides must be particularly aware of the risk of performing unsafe lifts or movements when placing the user on the patient/dinner chair. This applies whether the transfer is effected by means of a lift, slide board or support device.

MAINTENANCE AND CLEANING

The frame components and the seat and armrests may be cleaned with water and a household detergent. To thoroughly clean the seat: remove the seat, clean and make sure to fasten/press down the seat properly again.

Regularly check all bolts and nuts, and retighten if required.

⚠ Failure to inspect the bolts as directed may cause functional failures.

WARRANTY

The XXL-Rehab Patient/Dinner chair is covered by one year warranty from the invoice date for materials and workmanship.

Products which have been used for purposes other than the intended use; or which have been subject to negligence, abuse, improper storage or handling, improper operation, unauthorized modifications or damages beyond normal wear and tear as determined by XXL-Rehab, are not covered by this warranty.

If weight capacity is exceeded, the warranty will be void.

Any unauthorized repairs to product / part, as well as tampering with any components, will void the warranty.

XXL-Rehab Bedside Patient/Dinner chair is manufactured and tested in accordance with the international DS/EN 12182:2000 standard.

07/2014/JP

MANUFACTURER

XXL-Rehab
Fuglebækvej 1D
DK-2770 Kastrup
Denmark
Tlf.: +45 70 25 25 22
Fax: +45 70 25 25 44
www.xxl-rehab.com

Maximum weight capacity is	325 kg
Total length	80-101 cm
Total outer width on armrests	74 / 84 cm
Total height	79-94 cm
Seat width	61 / 71 cm
Seat depth	40-52 cm
Seat height	44-59 cm
Height of backrest	37 cm
Distance from floor to seat underside	33-48 cm
Distance between armrests	65 / 75 cm
Seat to armrest	28 cm
Weight	16.6 / 17.6 kg
Frame	Powder-lacquered chromium molybdenum steel
Seat / back / armrests	Polyurethane
ISO - classification	ISO 09 33 03

DK

UK

Strongly engineered and functionally designed exclusively for bariatric users!

XXL-Rehab Minimaxx

XXL-Rehab Minimaxx
Push Motor

XXL-Rehab Manual
Wheelchair - The Eclipse

XXL-Rehab Manual
Wheelchair - The Eclipse Tilt

XXL-Rehab
HD Electric Wheelchair

XXL-Rehab Wheelchair Cushions

XXL-Rehab Shower Bench

XXL-Rehab
Bedside commode

XXL-Rehab
Shower Commode

XXL-Rehab
Shower Commode Wheeler

XXL-Rehab
Shower Commode Rise N Recline

XXL-Rehab
Shower Commode Tilt

XXL-Rehab Barton Transferstol

XXL-Rehab Patient/Dinner chair

XXL-Rehab Rollator King

XXL-Rehab Rollator Support

XXL-Rehab Stand Tall

XXL-Rehab Bariatric walking frame

XXL-Rehab Canes

XXL-Rehab Multifunktional Scale

XXL-Rehab Lift and Slings

XXL-Rehab 400 bariatric hospital bed

XXL-Rehab Home Care Bed

XXL-Rehab Floor raiser

XXL-Rehab
Comfort Foam Mattress

XXL-Rehab
Low Air Loss mattress

XXL-Rehab
Alternating mattress

XXL-Rehab Bariatric suit

BRUGSANVISNING · USER MANUAL

mauryathome.dk

07-2014 XXL-Rehab Patientstol / XXL-Rehab Patient/Dinner chair

XXL REHAB®

XXL-Rehab · Fuglebækvej 1D · DK-2770 Kastrup · Denmark
Telephone: +45 7025 2522 · Fax: +45 7025 2544
info@xxl-rehab.com · www.xxl-rehab.com